



785
HELIOS IV™

Empowering
NANO with PICO

October 2021




Transforming Technology to Enrich Your Life

TABLE OF CONTENTS

- KEY BENEFITS
- CONFIGURATION
- NEW GUI
- WHY 785NM WAVELENGTH?
- REDEFINING CUTTING EDGE TECHNOLOGY
- TECHNICAL SPECIFICATIONS
- CONSTANT STABLE OUTPUT ENERGY & PULSE DURATION
- PROPRIETARY DOE FRACTIONAL TECHNOLOGY
- COMPARISON
- COMPETITION
- APPLICATIONS
- CLINICAL RESULTS
- 785nm vs 755nm
- PUBLICATION
- MARKETING SUPPORT



KEY BENEFITS






CONFIGURATION



System



Handpieces

*optional



NEW GRAPHIC USER INTERFACE





WHY **785nm** WAVELENGTH?

- Optimized 785nm wavelength for treating benign pigmented lesions and skin brightening
- High stability implemented by LASEROPTEK's proprietary picosecond technology enable faster and safer treatment
- Stronger melanin absorption than 1064nm
- Good for blue and green tattoo ink removal
- Minimized PIH

Faster Result

Safer for Darker Skin

Optimized for Pigmented Lesions






REDEFINING CUTTING EDGE TECHNOLOGY

Unique Laser Resonator Design

- LASEROPTEK's unique design of a solid-state laser resonator enables 785nm true pico pulse and 1064/532nm nano pulses in one platform.
- LO's micro cavity design enabled the high-power Ti:Sapphire laser.



Simplified H785 Laser Resonator Design Scheme




TECHNICAL SPECIFICATIONS

	Nanosecond		Picosecond
Wavelength (nm)	1064	532	785
Pulse Energy (mJ)	100 – 1400 300 – 2000 (RTP)	20 – 500	10 – 200
	300 – 3000 (FR)		
Pulse Duration	5 – 10 ns 300 µs (FR)		600 ps
Repetition Rate (Hz)	1 – 10		
Spot Size (mm)	5 x 5	4 x 4*	5 x 5
	1 – 7 (Zoom)		1 x 1 – 7 x 7
	2 – 10 (Zoom Colli)		7 x 7 (Colli)
Handpiece	1064 FX Zoom Zoom Collimator	Dia FX 785 S 785 Zoom 785 Collimator	
Display	10.4" 1024 x 768 TFT LCD		
Dimension (mm)	298 (W) x 819 (D) x 936 (H)		
Weight (kg)	80		

OPTIONAL				
Wavelength (nm)	532	785	595	660
Pulse Energy (mJ)	20 – 500	30 – 200	300	200
Repetition Rate (Hz)	1 – 10			
Spot Size (mm)	4 x 4	8 x 8		
Handpiece	532 FX	Dia FX 785	595 Dye	660 Dye




CONSTANT STABLE OUTPUT ENERGY & PULSE DURATION




785nm

■ H4 ■ 1-10 ■ 11-20 ■ 21-30



1064nm


■ H4 ■ 1-10 ■ 11-20 ■ 21-30




532nm

■ H4 ■ 1-10 ■ 11-20 ■ 21-30


Std Dev for Output Energy
< 2%



785nm



1064nm



532nm

Std Dev for Pulse Duration
< 1%

Above output energy values were measured using an energy meter and pulse duration graphs using an oscilloscope



PROPRIETARY DOE FRACTIONAL TECHNOLOGY

LASEROPTEK adopts DOE technology to its fractional handpieces. DOE (Diffractive Optical Element) is designed to split a single laser beam into a predefined number of beams. Each microbeam has uniformed energy intensity and penetration depth. The characteristic of DOE is that the quality of the output laser beam is independent of that of the input beam helps the output laser beams uniform and stable all the time.

Dia FX 785 is available in 5mm x 5mm and 8mm x 8mm*.

Twist & Treat




*optional

Dia FX 785 S

Dia FX 785*

1064 FX


532 FX*



Dia FX 785's 49 DOE micro beams



Dia FX 785 Penetration Depth



Lv. I

Lv. II


Lv. III

*Tested on crystal



THERMAL LENS COMPENSATING TECHNOLOGY

LASEROPTEK's proprietary thermal lens compensating technology prevent not the only thermal lens but also thermal birefringence to keep a very uniform beam quality and high output energy stability.






MULTIPLE LASER MODES

RTP (Real Twin Pulse) Mode


Core technology to have the same energy between two identical pulses



FR Mode

Genesis (Quasi Long) technique mechanism of action

Quasi- Long Pulse (300µs) Nd:YAG Laser Beams



** Stimulates skin cells regeneration*



WHAT'S NEW FROM HELIOS III (Nd:YAG)

	HELIOS III	HELIOS 785
Pulse Energy	Max 1300mJ Min 500mJ (RTP/FR)	Max 1400mJ Min 300mJ (RTP/FR)
Spot Size	1 – 8mm	1 – 10mm
Handpiece	1064/532 FX Zoom Collimator	1064 FX, 532 FX* Zoom Zoom Collimator
Display	10.4" 800 x 600 TFT LCD	10.4" 1024 x 768 TFT LCD
Facilities		Energy meter implemented Improved Zoom sensing Improved energy output efficiency



Energy meter implemented to monitor output energy status



Improved Zoom handpiece's sensing type to minimize sensing errors (Image sensing type)



Re-designed to improve the energy output efficiency



*optional



COMPARISON

Specifications	HELIOS 785	A	B
Wavelength	785nm	785nm	755nm
Pulse Energy	Max 200mJ	Max 100mJ	Max 200mJ
Pulse Duration	600ps	300ps	750ps
Peak Power	0.33GW	0.3GW	0.26GW
Repetition Rate	1 – 10Hz	1 – 10Hz	
Handpieces	Dia FX 785 S 785 Zoom 785 Collimator	785nm Ti:Sapphire	755 Zoom Fixed H/P





COMPETITION

Model	HELIOS 785	Revlite SI	PicoSure	StarWalker QX	Q Plus C Evo	Hollywood Spectra
Manufacturer	LO	Cynosure	Cynosure	Fotona	Quanta	Lutronic
Laser Type	Pico/QSNY	QSNY/ Dye	QS Alexandrite	QSNY	QSNY +Ruby	QSNY
Wavelengths (nm)	1064/532/785 *595/660	1064/532/ 585/650	532/755	1064/532 *585/660	1064/532/694	1064/532 *585/660
Pulse Energy (mJ)	1064: ~1400 532: ~500 FR: ~3000 1064RTP: ~2000 785: ~200	1064: ~1600	755: 200 532: 20	1064: ~1600 532: ~600 1064 long: ~15000 532 long: ~2000	1064: ~1500 532: ~500 694: ~1200	1064: ~ 1200/1400 532: ~ 400 Spectra: ~1500
Pulse Duration	785: 600ps 5~10nsec FR: 300µs RTP: 5 - 10nsec	5~20ns <7ns (585/650 dye)	532: < 600ps 755: 750ps	5~20 ns Long: 600~50,000 µs	6n dbl pulse; OP: 6ns + 150µs + 6ns PT: 300µs	5-10ns 1064 Q-3/Q-4: 10- 20ns Spectra: 190µs
Features	785 picosecond 200mJ 0.33GW	Gold standard, Soft beam	Focus Lens Array	Long mode (0.6~50msec)	Ruby 694nm	4 pulses Spectra mode (0.3ms)
Handpieces	1064/532 Zoom, Zoom Colli, 785 Zoom, Colli, 1064 FX, Dia FX 785 *532 FX, Dia FX 785	SI, 532Lite, Dye 585, Dye 650	755 Zoom, Fixed 532, Fixed 755	R28, R29, R58, FS20A, FS20B, FS20C, FS50B *R585, R650		Dual Focused Dots, Fractional (MLA), Zoom Colli, Zoom *585/650 Dye
Released Year	2021	2013	2012	2017	2012	2021



APPLICATIONS

- **Skin toning and skin brightening**


- Skin brightening and lightening
- Skin rejuvenation
- Genesis

- **Treating Benign pigmented lesion**

- Café-au-lait birthmarks
- Solar lentiginos, and senile lentiginos
- Becker's nevi
- Freckles
- Melanochiae
- nevus spilus, nevus of ota, and abnom
- PIH, etc.

- **Tattoo Removal:**

- black, blue, red, sky blue, green, red, and purple





CLINICAL RESULTS

785 Zoom for Epidermal Pigmentation
Dia FX 785 for 785 Pico Toning, 4-5 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Zoom for Epidermal Pigmentation
Dia FX 785 for 785 Pico Toning, 4-5 TX



BEFORE


AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

Skin whitening, Rejuvenation, and PIH
785 Dia FX, zoom, 150-200mJ, 5 TX



BEFORE

AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

Skin whitening, Rejuvenation, and PIH
785 Dia FX, zoom, 150-200mJ, 5 TX



Insert caption here

BEFORE

AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

Skin whitening, Rejuvenation, and PIH
785 Dia FX, zoom, 150-200mJ, 5 TX



Insert caption here

BEFORE

AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning, Age Spots
785 Dia FX, zoom, 100-200mJ, 3 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning, ABNOM
785 Dia FX, 150-200mJ, 5-7 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning, Solar lentigo, Melasma
785 Dia FX, 150-200mJ and 1064 collimator 500-600mJ
Combination treatment



BEFORE

AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning_785 Dia FX 150-180mJ, 5 TX
Seborrheic keratosis_785 Zoom, 180-190mJ, 2 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning_785 Dia FX 150-180mJ, 5 TX
Seborrheic keratosis_785 Zoom, 180-190mJ. 2 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning, Whitening, Blemish
785 Dia FX 150-180mJ, 5 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic




CLINICAL RESULTS

785 Toning_785 Dia FX 150-180mJ, 5 TX
Age spot_785 Zoom 180-190mJ, 1 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning_785 Dia FX 150-180mJ, 5 TX
Freckles_785 zoom 100-120mJ, 1 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

785 Toning, Age spots, Blemishes
785 Dia FX, zoom, 100-200mJ, 3 TX



BEFORE



AFTER

Photos Courtesy of Seok-bae Seo, M.D., SAS Dermatologic Clinic



CLINICAL RESULTS

1064nm Toning



1064nm Toning X 785nm



Photos Courtesy of Jong Gu Kim, M.D., Cheonan Oracle Dermatology Clinic



CLINICAL RESULTS

Congenital pigmented nevus

Treated with 755nm,
7 TX for 3 years



Treated with 785nm, 4 TX,
a photo taken in 3 days after treatment



*Photos Courtesy of Jong Gu Kim, M.D., Cheonan Oracle
Dermatology Clinic*



785nm vs 755nm




https://www.dropbox.com/s/kf8w57yz6kd0q8e/Skin%20Response_785%20vs%20755.mp4?dl=0

Photos Courtesy of Jong Gu Kim, M.D., Cheonan Oracle Dermatology Clinic



PUBLICATION



<https://doi.org/10.3390/app112210539>

Applied Sciences, 2021

Engineering of Ti:Sapphire Lasers for Dermatology and Aesthetic Medicine

Abstract: This review describes new engineering solutions for Ti:Sapphire lasers obtained at Laseroptek during the development of laser devices for dermatology and aesthetic medicine.

HELIOS-4, based on the Ti:Sapphire laser, produces 300mJ, 0.5ns pulses at 785nm for tattoo removal. The development of the last laser required the detailed study of a generation process and the investigation of the factors responsible for the synchronization of the generation in Ti:Sapphire lasers with short (several millimeters) cavities.

- Aleksandr Tarasov
- Hong Chu

Laseroptek Co., Ltd., Gyeonggi-do, Korea



PUBLICATION



<https://doi.org/10.1111/dth.15240>

- **Title:** Split-face comparative trial of 785-nm picosecond Nd:YAG laser and precision cryotherapy combination treatment for facial benign pigmented lesions
- Published on Dermatology Therapy, 2021
- **Authors:** Jae Wan Park, Hye Sung Han, Young Gue Koh, Suk Bae Seo, Gun-Ho Kim, Kyui Young Park
- **Key Points:**
 - The whole face of all patients was treated with a 600ps pulse using a 785nm Nd:YAG laser (HELIOS IV 785)
 - A total of 1000 pulses of laser therapy were administered using a single-parameter DOE with a fluence of 0.2J/cm², spot size of 5mm x 5mm and frequency of 10Hz
 - HELIOS IV 785
 - Clinical improvement of pigmented lesions and satisfaction of subjects were higher on the laser and cryotherapy combination treatment side than on the laser-only side.



MARKETING SUPPORT



· Brochure



· Banner



· Patient Leaflet

Quick Operation Guide HELIOS IV

EPIDERMAL PIGMENTED LESION

MODE	WAVELENGTH	SPOT SIZE	FLUENCE	REPETITION RATE	SPACING	NUM. OF PASSES	TIME	ANESTHESIA	RECOVERY TIME
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2

DERMAL PIGMENTED LESION

MODE	WAVELENGTH	SPOT SIZE	FLUENCE	REPETITION RATE	SPACING	NUM. OF PASSES	TIME	ANESTHESIA	RECOVERY TIME
HELIOS IV	785nm	1000um	11-15	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	11-15	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	11-15	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	11-15	10-15	1-2	1-2	1-2	None	1-2

TATTOO REMOVAL

MODE	WAVELENGTH	SPOT SIZE	FLUENCE	REPETITION RATE	SPACING	NUM. OF PASSES	TIME	ANESTHESIA	RECOVERY TIME
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2
HELIOS IV	785nm	1000um	14-20	10-15	1-2	1-2	1-2	None	1-2

LASEROPTIK

· Quick Guide



MARKETING SUPPORT



· H785 Promo Video



THANK YOU

